### SHORT REPORT Open Access

## CrossMark

# Extraction of total nucleic acids from bacterial isolates using the bioMérieux NucliSENS easyMAG total nucleic acid extractor

Eleanor A. Powell and Joel E. Mortensen\*

#### Abstract

The BioMerieux NucliSENS easyMAG total nucleic acid extractor was evaluated for use on bacterial isolates in the clinical microbiology laboratory. Forty eight isolates were extracted, yielding quantifiable amounts of DNA for all isolates. The easyMAG is appropriate for DNA extraction from bacterial isolates and will be incorporated in the clinical laboratory.

**Keywords:** DNA extraction, 16S rRNA sequencing, Bacterial sequencing

The number of molecular methods and their importance in the clinical microbiology laboratory is rapidly increasing. Nucleic acid extraction from bacterial isolates is an important first step in several molecular test methodologies, including 16S rRNA sequencing for bacterial identification, whole genome sequencing, and other genotypic methods. While some clinical microbiology testing platforms have integrated nucleic acid extraction, others require separate nucleic acid extraction before testing may begin. Laboratory quality control has recently come under renewed scrutiny with the implementation of individual quality control plan (IQCP) requirements and as part of this, it is critically important to identify steps in a testing process that may result in testing failure or testing delays.

Inadequate nucleic acid extraction may result in an insufficient amount of nucleic acids, which can result in wasted reagents and delayed results. Manual extraction methods have traditionally been used, but these methods require extensive hands-on time and frequently have many steps, which increases the opportunity for method error or sample contamination [1]. The NucliSENS

nucleic acid extraction from various sources, including blood, stool, and urine. The easyMAG extracts nucleic acids using magnetic silica beads. EasyMAG extraction of clinical samples such as respiratory samples, stool, and blood followed by molecular detection of viral and bacterial pathogens has been reported [2-4]. While the easyMAG has been used on a variety of clinical specimen types, the extractor does not include pre-set protocols for bacterial extraction. Although a small number of published studies have utilized the easyMag extractor to extract nucleic acids from bacterial isolates, these have typically been limited to a small range of species [5]. For the easyMAg to be utilized as the primary bacterial nucleic acid method in a clinical microbiology laboratory, the extractor must be capable of extracting nucleic acids from a wide variety of bacterial species. To determine the suitability of the easyMAG for routine bacterial nucleic acid extraction, the current study sought to examine the ability of the easyMag extractor to extract nucleic acids from a wide variety of bacterial isolates.

easyMAG total nucleic acid extractor (bioMérieux,

Marcy l'Etoile, France) is a platform capable of total

Forty-eight isolates from the American Type Culture Collection (ATCC) were selected to represent the wide variety of isolates typically encountered in a clinical microbiology laboratory (Table 1). Organisms were

<sup>\*</sup>Correspondence: joel.mortensen@cchmc.org Diagnostic Infectious Disease Testing Laboratory, Cincinnati Children's Hospital Laboratory, 3333 Burnet Avenue, MLC 1010, Cincinnati, OH 45229, USA



Table 1 Bacterial species included in validation

Gram-positive cocci	Gram-positive bacilli	Enterobacteriaceae	Non-Enterobacteriaceae gram-negative organisms	
Staphylococcus epidermidis Streptomyces albus		Escherichia coli (3)	Campylobacter jejuni	
Staphylococcus aureus	Exiguobacterium mexicanum	Providencia stuartii	Campylobacter coli	
Enterococcus faecium	Clostridium septicum	Salmonella enterica	Pseudomonas aeruginosa (2)	
Streptococcus pyogenes	Actinomyces viscosus	Shigella sonnei	Prevotella melaninogenica	
Staphylococcus saprophyticus	Propionibacterium acnes	Plesiomonas shigelloides	Stenotrophomonas maltophilia	
Enterococcus faecalis	Corynebacterium striatum	Enterobacter cloacae	Alcaligenes faecalis	
Streptococcus pneumoniae	Clostridium perfringens	Vibrio vulnificus	Acinetobacter baumannii	
Streptococcus mitis	Actinomyces pyogenes	Klebsiella oxytoca	Haemophilus influenzae	
Streptococcus agalactiae	Streptomyces griseus	Citrobacter freundii	Burkholderia cepacia	
	Listeria monocytogenes	Proteus mirabilis (2)	Neisserisa meningitidis	
		Klebsiella pneumoniae	Pseudomonas putrefaciens	
		Serratia marcescens	Pasteurella multocida	
			Bordetella bronchiseptica	
			Moraxella catarrhalis	

Table 2 Concentration and A260/A280 of extracted total nucleic acid

	DNA concentration (ng/μL)		A <sub>260</sub> /A <sub>280</sub>	
	Median	Range	Median	Range
Gram-positive cocci	8.5	7.4–14.5	1.35	1.05-1.70
Gram-positive bacilli	9.1	8.0-37.8	1.47	1.09-1.80
Enterobacteriaceae	15.9	12.1-20.4	1.63	1.27-1.94
Non-Enterobacteriaceae	17.4	7.8-22.0	1.67	1.02-1.82

assigned to one of four groups: Gram-positive cocci, Gram-positive bacilli, Gram-negative bacilli in the family *Enterobacteriaceae*, and Gram-negative bacilli other than *Enterobacteriaceae*. A 1.0 McFarland suspension of each isolate in 0.45 % saline was made using the bioMérieux DensiCheck. 100 microlitre of each isolate was placed in the vessel portion of the easyMAG disposables and inserted in the instrument. DNA extraction was set up using the following parameters: program-specific 1.0.1, volume (input)-0.1 mL, elution volume-110  $\mu$ L; lysis- on board. Ten microlitre of the resulting nucleic acid eluent was quantified using the Eppendorf Biophotometer (Eppendorf, Hauppauge, New York, USA).

For each of the four groups, the two isolates with the lowest nucleic acid concentrations were submitted for 16S rRNA sequencing. When these two isolates were from the same genus, the isolate with the next lowest nucleic acid concentration was submitted, ensuring more diverse representation of each group. Sequencing results were compared to the ATCC identification of the organism to ensure correct identification. Four additional

isolates yielding higher extracted DNA concentrations  $(7.8-21.6 \text{ ng/}\mu\text{L})$  were also sequenced.

Extraction with the easyMAG resulted in detectable amounts of nucleic acids for all isolates (Table 2). Grampositive organisms generally resulted in lower nucleic acid concentrations and A<sub>260</sub>/A<sub>280</sub> ratios than Gramnegative organisms.  $A_{260}/A_{280}$  is a measure of nucleic acid purity [6].  $A_{260}/A_{280}$  of 1.8 represents pure DNA, while an A260/A280 of 2.0 represents pure RNA. Lower values are typically the result of protein contamination. For those isolates with the lowest concentrations of nucleic acids in each group, 16S rRNA sequencing was successful. It is expected that isolates yielding higher concentrations of DNA would also sequence, as low DNA quantity is frequently cited as the reason isolates cannot be sequenced. To confirm this, four isolates with higher extracted DNA concentrations were also sequenced, and all yielded identifications. Low A<sub>260</sub>/A<sub>280</sub> also did not affect downstream 16S rRNA testing, as isolates with  $\rm A_{260}/A_{280}$  as low as 1.05 were sequenced successfully. Low ratios are likely the result of protein in the DNA extract, which does not prevent downstream testing. It is expected that nucleic acids from easyMAG extractions will be suitable for other PCR-based methodologies.

In-house nucleic acid extraction can decrease turnaround-time of downstream testing and minimizes the risks associated with sending potential pathogens to other laboratories. Coupling laboratory-performed nucleic acid extraction with nucleic acid quantification allows laboratories to perform additional manipulation on isolates that initially do not yield adequate amounts of nucleic acid. Identification of these isolates before downstream testing can prevent unsuccessful testing, save reagents and testing time, and decrease overall time to test result. While this was not required in the current study, lysis can be performed manually rather than on the instrument to extend lysis time and additional physical lysis steps such as grinding can also be added. Additional studies are necessary to develop protocols for nucleic acid extraction from acid-fast bacilli and fungi. The efficacy, flexibility, and ease-of-use of the bioMérieux NucliSENS easyMAG total nucleic acid extractor make it suitable for use in the routine clinical microbiology laboratory.

#### Authors' contributions

EAP and JEM made substantial contributions to the study design, acquisition and analysis of the data. Both authors read and approved the final manuscript.

#### Competing interests

The authors declare that they have no competing interests.

#### **Funding**

No outside funding sources.

Received: 17 May 2016 Accepted: 13 September 2016

Published online: 22 September 2016

#### References

- Loens K, Bergs K, Ursi D, Goossens H, Leven M. Evaluation of NucliSens easyMAG for automated nucleic acid extraction from various clinical specimens. J Clin Microbiol. 2007;45(2):421–5.
- Esona MD, McDonald S, Kamili S, Kerin T, Gautam R, Bowen MD. Comparative evaluation of commercially available manual and automated nucleic acid extraction methods for rotavirus RNA detection in stools. J Virol Methods. 2013;194(1–2):242–9.
- Kalina WV, Douglas CE, Coyne SR, Minogue TD. Comparative assessment of automated nucleic acid sample extraction equipment for biothreat agents. J Clin Microbiol. 2014;52(4):1232–4.
- Manji R, Zheng X, Patel A, Kowerska M, Vossinas M, Drain A, Todd KM, Lenny N, DeVincenzo JP, Ginocchio CC. Multi-center evaluation of the adenovirus R-gene US assay for the detection of adenovirus in respiratory samples. J Clin Virol. 2014;60(2):90–5.
- Monteiro J, Widen RH, Pignatari AC, Kubasek C, Silbert S. Rapid detection of carbapenemase genes by multiplex real-time PCR. J Antimicrob Chemother. 2012;67(4):906–9.
- Manchester KL. Value of A260/A280 ratios for measurement of purity of nucleic acids. Biotechniques. 1995;19(2):208–10.

# Submit your next manuscript to BioMed Central and we will help you at every step:

- We accept pre-submission inquiries
- Our selector tool helps you to find the most relevant journal
- We provide round the clock customer support
- Convenient online submission
- · Thorough peer review
- Inclusion in PubMed and all major indexing services
- Maximum visibility for your research

Submit your manuscript at www.biomedcentral.com/submit

